

High Performance Solutions

Forefront Medical Technology's enteral feeding tubes are a proven choice for reliable patient feeding. They are constructed from high performance, medical grade, polyurethane tubing, and manufactured with a Hydromer® water-activated coating that eliminates the need for additional lubrication. The precision extruded, thin wall tubing design provides a higher delivery flow rate of nutritional formula to the patient. Custom connectors are available to meet OEM specific requirements.

The device package includes a stainless steel stylet with a unique spring cushion tip, that facilitates easy insertion into the patient. The dual port connector provides for both simultaneous feeding and decompression. Special connector designs are available for unique design protocols and applications. Forefront's feeding tubes are designed to resist kinking or collapse, during aspiration. The tubes are radiopaque and include reference markings.



Specialty Designs for Unique Applications

Forefront Medical's standard products include 8 to 12 FR tube diameters, and are available in two types: Naso-Gastric (NG) and Jejunal (J) feeding tubes. Typical specifications for our standard product range are summarized in the table below. Jejunal feeding tube catheters are available for both guide wire and tether placement techniques. Special designs, sizes, lengths, weight sets or other unique requirements are available upon request.

Forefront Medical Technology is the integrated global contract manufacturer for disposable, diagnostic and drug delivery devices. Please contact us to confidentially discuss your application at (860) 255.7610 (USA) or (86) 21 6062 7177 (Shanghai), or visit our website www.forefrontmedicaltechnology.com.

Product #	Size	Description	Overall Length (in.)	OD (in.)	ID (in.)	Insert
Naso-Gastric Catheters (NG Tubes)						
NG6896	8 FR	Enteral Feeding Tube	45	0.115	0.083	Stylet Wire
NG6897	10 FR	Enteral Feeding Tube	45	0.135	0.094	Stylet Wire
NG6898	12 FR	Enteral Feeding Tube	45	0.157	0.109	Stylet Wire
Jejunal Feeding / Gastric Decompression Catheters (J Tubes)						
JT7099	9 FR	Suture Tip J Tube (Pull)	38.37	0.118	0.080	Stylet Wire
JT7098	9 FR	Wire Guided J Tube (Push)	37.13	0.118	0.080	Guidewire
JT7101	12 FR	Suture Tip J Tube (Pull)	37.93	0.157	0.109	Stylet Wire
JT7102	12 FR	Wire Guided J Tube (Push)	37.13	0.157	0.109	Guidewire

